ES

Nº registro: 1648 RD

Reg. ICA 10194 BV

Reg Sanitario:

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EN

Giardia
Diagnostic Kit
Uranotest®
For veterinary use only

Technical basis

The UranoTest® Giardia diagnostic kit is based on the immuno chromatographic technique and is designed for the qualitative detection of antigens of Giardia in feline and canine faeces, with a high sensitivity.

The test consists of several overlapped membranes. On one of the membranes there is a test line (T) and a control line (C). The lines are not visible before the test begins. Once the test begins, the sample migrates through the membrane until it reaches the control line. If the result is negative, only the control line will appear. If the result is positive, both the test and control lines will appear. The migration is driven by capillary action through the membrane.

Materials supplied

- 1 - Test devices individually packaged in aluminium pouch.
- 2 - Tubes with buffer solution for sample dilution.
- 3 - Swabs for sample collection.
- 4 - Disposable pipettes.
- 5 - Instructions for use.

Precautions

- 1 - For veterinary use only
- 2 - Wear disposable gloves when handling the samples. All samples should be treated as potentially infectious. Wash and disinfect hands after handling.
- 3 - All samples must be handled as potentially infectious.
- 4 - Avoid aerosol formation when dispensing the sample.
- 5 - To obtain good results, it is important to add the correct sample volume.
- 6 - Do not use the test if the envelope is damaged or broken.
- 7 - Do not reuse.
- 8 - All reagents must be at room temperature before performing the test.
- 9 - Do not use reagents after the expiry date.
- 10 - The quality of each component of the kit has been individually assessed for each batch. Do not mix components or reagents from kits with different batch numbers.

Preservation and stability

The kit must be stored at a temperature between 2 and 30° C. Under these conditions, we can guarantee the stability until the expiry date printed on the box and on the individual pouch.


The kit has been developed to be stored at room temperature. Although it also can be stored in the refrigerator, we recommend storing it at room temperature to avoid the need to wait for reagents to reach the room temperature.

DO NOT FREEZE. Do not expose to direct sunlight.

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
Sample collection

- 1 - Dog and cat faeces can be used.
- 2 - The sample may be collected directly from the rectum using the provided swab. Add sufficient amount of sample to carry out properly the test is shown in the picture.
- 3 - Insert the swab into the tube containing the diluent and press it against the tube walls in order to release the sample that could contain the virus.
- 4 - Shake the tube in order to ensure a good homogenization.

**Instructions for use**

- 1 - Remove the test device from the protective pouch and place it on a flat and dry surface.
- 2 - Using the pipette provided, transfer 4 drops of the recently prepared sample to the test area. Note that the sample must contain particles, wait for 1 minute until they sediment, use the supernatant.
- 3 - When the test begins you will observe migration of the sample through the result window. If migration has not begun 1 minute after sample addition, add one more drop of the diluted sample.
- 4 - Read the results within 5-10 minutes. Coloured lines appeared after 20 minutes have not diagnostic value and should be ignored.

**Interpreting results**

- 1 - Negative result
There is only a single purple line on the C area. This line should always appear.
- 2 - Positive result
Two lines appear on the result window. Whichever line appears first, the result is considered positive.
- 3 - Invalid result
The test is invalid if no coloured line appears in the Control area (C) even if a coloured line appears in the Test area (T).

Limitations of the technique

Even though the UranoTest® Giardia Vet kit shows high sensitivity and specificity, cannot be excluded a low incidence of false positive or negative results. As any other laboratory procedure, the definitive clinical diagnosis cannot be based only on the test result. It must be based on an ensemble of clinical and laboratory procedures. If there is any doubt, repeat the test and/or contrast with other diagnostic methods.

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